

2023 CASE STUDY

Precision Ethanol Recovery for Pharmaceutical Innovation



Industry: Pharmaceutical & Biotech

Unit Installed: Custom Fractional Distillation Column System

The Challenge

A pharmaceutical innovator needed a high-throughput ethanol recovery solution to support its biomedical production processes. Their feedstock consisted of a 90% ethanol and 10% water mixture, requiring precise separation to meet purity standards.

Key challenges included:

- Maintaining pharma-grade purity with 95% ethanol output.
- Processing over 100 gallons per hour continuously.
- Ensuring compliance with C1D2 safety standards.
- Integrating solvent management tools for tracking and weighing.

Traditional systems lacked the automation, throughput, and traceability the client required for regulated pharmaceutical operations.

maratek.com est. 1967

The Solution

Maratek engineered and installed a custom fractional distillation column system, featuring:

- 25-ft tall packed column with 316 stainless steel construction for sanitary, pharma-grade processing.
- Automated PLC controls (Allen Bradley or Siemens) with remote access and real-time monitoring.
- Solvent management system with integrated weigh scale and tracking for compliance and inventory control.
- High-capacity throughput of ~5,000L/day (24-hour operation).
- C1D2 certification for safe operation in hazardous environments.
- Onsite installation and training with two Maratek technicians for up to 7 days, ensuring smooth startup and operator confidence.

The Result

The client achieved:

- Consistent ethanol purity at 95%, meeting stringent pharmaceutical standards.
- · Operational efficiency with automated, high-volume recovery.
- · Improved traceability through integrated solvent tracking and weighing.
- Regulatory confidence with C1D2-certified equipment and full documentation.
- Seamless onboarding thanks to Maratek's hands-on training and support.

Why Maratek?

With over 50 years of experience and 1,500+ systems installed, Maratek is trusted for its custom engineering, automation expertise, and industry-leading safety standards. For this client, Maratek delivered a solution that aligned with pharmaceutical-grade requirements while enabling scalable, compliant operations.

Call: 905.857.2738 | Toll Free: 1.800.667.6272 | Email: sales@maratek.com

maratek.com est. 1967